

**CAN'T
HEAR?**

Our listening products will help you to hear better, whether you have hearing aids or not. They amplify sound and block out background noise.

Personal listeners

These are fantastic portable products that help you to hear more clearly, whether you have hearing aids or not. You can carry them around with you, and point them in the direction you want to listen. They amplify the sound and transfer it directly to your hearing aids or ears, via a variety of options:

- Headphones, earbuds and stethosets (if you don't have hearing aids).
- Ear hooks and neckloops (if you have hearing aids with a T setting).

Benefits of personal listeners

- They enhance the volume and clarity of conversations and other sounds.
- You can use them in most situations including chatting with friends, watching TV and listening in noisy environments.
- On some models you can adjust the high and low tones (bass and treble) to suit your particular hearing loss.
- Listeners with a T setting allow you to benefit from loop systems.

Room loops

Loops are generally for people who wear hearing aids with a T setting. They normally consist of an amplifier and a cable that goes around the perimeter of your room and plugs back into the amplifier. The amplifier then plugs into a sound source like a TV and

you can also attach a small microphone for conversations. The electrical current in the loop produces a magnetic field that sends sound directly to your hearing aids when they are switched to T.

How room loops can help you:

- They only amplify the sound you want to hear, whether that's the TV or a friend speaking into the microphone.
- There is no limit to the number of people that can use a room loop at one time.
- Organisations can use commercial loop systems so customers and staff can hear more clearly.

Elite listeners

Elite listeners have all the benefits of a personal listener, plus a transmitter.

Additional benefits


- The transmitter means you can hear someone from a long distance (at least 25m). The person speaking talks into the transmitter and the sound is picked up by the receiver.
- Some of our listeners allow you to change the direction from which sounds are picked up by the microphone:

Omni-directional - the microphone picks up sounds from all directions. This is useful in situations like meetings or dinner parties.


Directional - you can point the receiver in the direction you want to hear, and the sounds from other directions are made quieter.

If you wear a pacemaker, we advise you to speak with your GP and your pacemaker supplier before using a neckloop.


- A** Headphones
- B** Earbuds
- C** Stethoset
- D** Neckloop
- E** Ear hooks




Examples of an elite listener in use




In a café or restaurant




Transmitting from a distance




At a meeting



In the evening



Listening to a speaker



Out shopping



L335A/L389/L329

Action on Hearing Loss Sonido digital listener

The Sonido is excellent value as it gives you digital sound similar to that used in digital hearing aids, but at a fraction of the price of other digital listeners. Sonido isn't just digitally controlled – it has full digital processing, which improves sound quality and reduces background noise and feedback.

Conversation mode

The built-in microphone is automatically switched off when you connect the Sonido to your TV or stereo. So when you watch TV you can't hear the person sat next to you - but if you press the Mic button you will be able to hear the TV and the conversation simultaneously, while the button is pressed. When you release the Mic button you will return to the plug-in sound.

- Excellent for people who have trouble holding small things – it has an easy grip rubberised coating and large, raised, soft-touch buttons
- External microphone (included) allows you to amplify sounds directly from other sources, such as TV, radio or MP3 player
- Volume and tone control – useful for tailoring the bass and treble to your hearing loss

- 150 hours' battery life – almost double that of most other listeners
- Hinged battery cover for easy access
- Battery status light
- Clear LEDs show volume and tone
- To protect your hearing, Sonido automatically switches to a lower volume setting when turned on
- You can set the volume at different levels for each ear via a button in the battery compartment
- T setting that allows people who don't wear hearing aids to benefit from induction loops or test they are working
- Size: 137 x 45 x 25mm

First choose your Sonido colour and then choose how you would like to listen.

Options available

L376Blue

£56.68 exc VAT

£68.02 inc VAT

L376Silver

£56.68 exc VAT

£68.02 inc VAT

2 x AA batteries included with both options

You need one of the following products to go with your Sonido.

L333A Headphones

£14.78 not VAT exempt

L334A Stethoset

£14.48 not VAT exempt

L335A Neckloop

£12.50 exc VAT

£15.00 inc VAT

L329 Dual earhooks

£18.75 exc VAT

£22.50 inc VAT

L389 Single earhook

£18.40 exc VAT

£22.08 inc VAT

Accessory

L333B TV kit – 5m lead

£23.94 not VAT exempt

TV kit includes phono lead, scart adaptor and jack adaptor

It is excellent for the price. All the controls are straightforward and the sound quality is very good.
G F Cutter, Powys



Bellman Mino

This pocket-sized personal listener is very light and has digital sound processing, which is normally reserved for elite listeners. Digital sound means you get better sound quality and reduced background noise. The Mino is particularly good at filtering out continual sounds, such as the hum of a fan, while enhancing the sound of speech. The directional microphone amplifies sounds from the front while suppressing noise from the side and rear. You can also switch to the omni-directional microphone mode, which enables you to hear clearly from all directions: a useful feature in meetings or at dinner parties. The adjustable tone control means you can adjust the bass and treble to suit your particular hearing loss.

- Adjustable tone and volume control
- Reduces background noise and enhances speech
- Directional and omni-directional microphone settings
- Single use buttons make changing functions simple
- Up to 18 hours' operation from 2.5 hour charge
- Balance control, so you can adjust the volume for each ear (except with neckloop)
- Size: 99 x 48 x 22mm

Options available

L394 Mino with neckloop

£165.00 exc VAT

£198.00 inc VAT

L395 Mino with headphones

£165.00 exc VAT

£198.00 inc VAT

L396 Mino with earphones

£165.00 exc VAT

£198.00 inc VAT

Accessories

L333B TV kit - 5m lead

£23.94 not VAT exempt

TV kits includes phono lead, scart adaptor and jack adaptor

L394A Microset - earphone and lapel microphone

£20.94 not VAT exempt

This enables you to pop the Mino in your pocket and clip the microphone to your clothing



Start to sign! package

Our multimedia package Start to sign! is an essential support aid to class-led courses in British Sign Language (BSL) and an ideal introduction to learning. Now in its fourth edition, the book is bigger and brighter, with brand new colour photographs of more than 800 signs.

PLUS! Buy the book and get FREE access to our Start to sign! website. You'll be able to view more than 1,000 video clips - and test yourself with our online quizzes.

RN/2535/09/07

Start to sign! package

£19.99 not VAT chargeable

Works best on PCs with Internet Explorer or Apple Macs with Firefox



The Pocket Talker

This is a powerful portable listening device with a plug-in microphone. It makes listening to conversations easier whether in the car, at a restaurant, in small groups or at meetings. It's also very useful for listening to the TV, radio or hi-fi, and fits neatly into a top pocket.

- Tone control
- Five-year warranty on the amplifier only
- Choice of neckloop, over-ear earpiece, single earphone or headphone connection
- Noise-cancelling function
- Handy belt clip
- Battery life of up to 100 hours
- Size: 95 x 65 x 25mm



Options available

AP44 Pocket Talker with headphones/earphone/microphone extension lead

£66.50 exc VAT
£79.80 inc VAT

AP55 Pocket Talker with neckloop

£81.95 exc VAT
£98.34 inc VAT

AP40 Pocket Talker with over-ear earpiece

£66.50 exc VAT
£79.80 inc VAT

Accessories

AP43E Optional/spare headphones

£12.52 not VAT exempt

AC11 Optional/spare neckloop

£32.00 exc VAT
£38.40 inc VAT

AP43B Optional microphone extension cord (2.5m) with TV clips

£10.43 not VAT exempt

2 x AAA batteries included



Crescendo 50

The Crescendo is a compact, lightweight listener with a built-in microphone and removable belt clip so you can carry it easily. It has an easy-to-use variable tone and volume control feature. This allows you to manage the sound to suit your needs in different situations, such as listening to conversations, the television, hi-fi or radio. There's even a loop listening facility for people without hearing aids to use loop systems.

- Rotary volume and tone control
- Supplied with 4m tie-clip microphone and lead
- Battery life of up to 75 hours
- Size: 95 x 62 x 25mm

Options available

L322 Crescendo 50/1 with stethoset
£52.50 exc VAT **£63.00** inc VAT

L323 Crescendo 50/5 with neckloop
£52.50 exc VAT **£63.00** inc VAT

2 x AA batteries included

Accessories

X613 Replacement stethoset
£15.00 not VAT exempt

L335A Replacement neckloop
£12.50 exc VAT **£15.00** inc VAT

X672 Replacement tie-clip microphone (4m lead)
£11.88 not VAT exempt

L322A Scart lead (5m)
£11.34 not VAT exempt



MiniTech+

This easy-to-use portable listener is small and lightweight, so you can use it in the home or out and about. You can enhance a wide range of sounds easily, quickly and accurately using the volume control. The tone control lets you adjust the bass and treble to suit your hearing. It can also be used with a scart lead to listen to the TV. Comes with stethoset, neckloop, earbuds, tie-clip microphone and case.

- Push button volume control
- Tone control
- Ideal for people who find holding things difficult
- On/off and battery status indication
- Battery life of up to 90 hours with stethoclips or headphones (40 hours with a neckloop)
- Size: 95 x 55 x 22mm

L400 Minitech+ with stethoset, neckloop and earbuds

£58.00 exc VAT **£69.60** inc VAT
 2 x AA batteries included

Accessories

L292A Replacement stethoset
£18.00 not VAT exempt

L291A Replacement neckloop
£15.00 exc VAT **£18.00** inc VAT

L291B Replacement extension lead (3m) for a microphone
£18.00 not VAT exempt

L291C Replacement stubby microphone
£18.00 not VAT exempt

L291D TV kit (scart, phono and jack)
£18.00 not VAT exempt



Phonak EasyLink+

EasyLink is the perfect introduction to wireless personal listeners. It only has one button so it couldn't be easier to use. This is a good choice if you want the ideal combination of simplicity and the fantastic sound quality you would expect from a leading hearing aid manufacturer like Phonak.

- Omni-directional microphone helps you hear sounds from all directions
- One button operation – it simply turns on and off
- LCD screen showing FM channel and battery status
- Charging time of just two hours, with up to 12 hours' consistent use from one charge
- You can plug into your TV, MP3 player, stereo or PC for a clearer sound
- Operating range of up to 30 metres
- Can be connected to analogue telephones (see accessory page 41) helping you hear calls more clearly
- Size: 105 x 35 x 18mm (Receiver) 66 x 29 x 12mm (Transmitter)

L369 EasyLink+ and MyLink package

£500.00 exc VAT

£600.00 inc VAT



Phonak ZoomLink+

The ZoomLink has all the features of the EasyLink plus the following additional benefits. The ZoomLink has three microphone settings, which mean it's flexible enough to help you hear in a wide variety of situations:

Super zoom – allows you to focus the sounds in a narrow position in front of the transmitter – ideal for one-to-one situations.

Zoom – this mode allows you to pick up sounds in a specific direction – ideal for listening in a small group.

Omni – allows you to pick up sounds all around (360°) the transmitter – ideal for large gatherings.

- Three microphone settings for hearing in different situations
- Size: 105 x 35 x 18mm (Receiver) 66 x 29 x 12mm (Transmitter)

L368 ZoomLink+ and MyLink package

£550.00 exc VAT

£660.00 inc VAT





Phonak SmartLink+

SmartLink+ is the most sophisticated listener in the Phonak range. It is Bluetooth compatible, so you can use any Bluetooth mobile phone or Bluetooth audio device (some laptops and MP3 players) in a range of up to 10m. If you are listening to music on your MP3 player and your phone rings, the SmartLink+ will instantly pause the music and send the phone's ringing signal to the MyLink receiver. You simply

accept the incoming call by pressing the green telephone button on the SmartLink+ and talking directly into the transmitter. To end the call press the red button, and the MP3 player will automatically resume playing the music.

React to different situations at the press of a button

At the press of a button you can switch between three microphone options. The three microphone settings enable you to adjust how much of your environment you want to hear:

Super zoom - Allows you to focus on sounds within a narrow area in front of the transmitter. This setting is particularly good at enhancing speech and blocking out unwanted sounds in extremely noisy situations.

Zoom - This mode is similar to 'super zoom' but the microphone will pick up sounds from a wider area. This will filter out unwanted sounds and improve the speech signal.

Omni - Allows you to pick up on sounds all around (360°) the transmitter. This is particularly useful in meetings: you can put the SmartLink+ in the centre of the table and hear what

everyone is saying, rather than just one or two people.

- Bluetooth compatible
- Can be used with your telephone (see accessory below), Bluetooth mobile, laptop and MP3 player
- Rechargeable battery with up to 12 hours' use from two hours' charge
- Belt clip
- Keypad lock to prevent controls being pressed by accident
- Operating range between receiver and transmitter of up to 30m
- Operating range of the Bluetooth function is 10m

L370 SmartLink+ and MyLink package

£600.00 exc VAT

£720.00 inc VAT

Accessories

L401 Telephone adaptor

£24.97 not VAT exempt

If you pay by credit card you can spread the cost of your order over four months, if the face value is more than £100. Just tick the relevant box on the order form and your card will be debited each month in four equal amounts.





The Conversor Pro

The Conversor Pro is a discreet radio microphone with a pendant neckloop receiver that transmits sound to your hearing aids. The powerful microphone can be handheld or placed on a surface near the sound source. It gives you the ability to 'focus' and choose what you hear, whether in a one-to-one situation, a group discussion or a lecture hall.

The Conversor Pro reduces confusing background noise and lets you take part comfortably in most everyday conversations by improving sound intelligibility. The package also includes a neck cord which you can plug into the microphone to give you hands-free operation. This is ideal for teachers or lecturers to use during lessons, and is also great for someone with you to use while driving or shopping.

- Neckloop receiver has volume control
- Microphone boost button
- Headphone socket on receiver
- Rechargeable battery with up to eight hours' use from three hours' charge
- TV connection kit supplied (jack-to-jack lead, phono adapter lead and scart plug)
- Operating range of up to 25m
- Size: 100 x 40 x 35mm (transmitter) 65 x 40 x 25mm (receiver)

L343 Conversor Pro
£349.00 exc VAT
£418.80 inc VAT

Accessories

L333A Optional headphones
£14.78 not VAT exempt

L343A Replacement transmitter*
£165.00 exc VAT
£198.00 inc VAT

L343B Additional receiver*
£175.00 exc VAT
£210.00 inc VAT

L343C Replacement battery
£13.20 not VAT exempt

* Please note we will need to know what channel you are using on your existing system – A, B or C.

If you pay by credit card you can spread the cost of your order over four months, if the value is more than £100. Just tick the relevant box on the order form and your card will be debited each month in four equal amounts.





Bellman Domino Classic

Domino Classic uses digital state-of-the-art sound processing to enhance speech while reducing distracting background noise. You can also adjust the tone of people's speech so that the sound suits your particular type of hearing loss.

It's simple to use - all the common features can be accessed with a single push of a button.

You can hear sounds through the transmitter microphone and the receiver microphone and balance the sound between them. This is useful in situations where you want to clearly hear two people at once. You can easily control all the functions of the two units through the receiver, so there is no need to reach for a distant transmitter when you want to make a volume or tone adjustment. It comes with a durable travel case, together with an international USB-powered charger kit.

- Volume and tone control
- Omni-directional microphone in both the receiver and transmitter
- Domino uses digital wideband technology to stream audio in full stereo. The result is crystal-clear stereo sound with outstanding clarity and definition

- State-of-the-art digital encryption technology makes it impossible to eavesdrop
- Works on a frequency that can be used anywhere in the world
- No handling noise, due to a special internally suspended microphone
- You can listen to the TV or a conversation, or both at the same time
- You can use up to eight receivers with each transmitter
- Built-in rechargeable battery with 8-12 hours' use from 2.5 hours charge
- Operating range of up to 25m inside and 150m outside
- Size: 99 x 48 x 22mm

L393 Domino Classic with neckloop and earbuds

£399.00 exc VAT

£478.80 inc VAT

Just what I wanted. I can listen to the TV and hear the phone whilst on loop.

Paul Upton, Exmouth



Bellman Domino Pro

This high-spec elite listener has all the features of the Domino Classic plus the following additional benefits:

- The microphones in the receiver and transmitter give greater flexibility as you get a directional option as well as the omni-directional option of the Classic
- The balance control function is useful if you have a different hearing level in each ear
- You have a choice of mono or stereo - mono is sometimes preferred by people with one hearing aid, so that they hear all the sound in that ear
- A tie-clip microphone can be attached to a transmitter, a discreet method for presentations
- A connection lead plugs into your TV, stereo or MP3 player, giving you the best possible sound with no external noise whatsoever

L383 Domino Pro with neckloop and earbuds

£625.00 exc VAT

£750.00 inc VAT



L320 only

2 year warranty

Contego digital listener

This stylish, lightweight and easy-to-use wireless listener is truly multi-functional. You can use it for TV listening, conversations with friends, in a meeting or lecture, in the car or even while you're shopping. Whoever is speaking wears the microphone, which transmits the sound to a receiver worn by the user. You can be up to 25m away and still pick up the speaker's voice. Hearing aid users with a T setting can use it with the neckloop option. If you don't wear hearing aids, you could choose our other option which includes headphones and earphones.

The receiver has its own built-in microphone, which allows you to not only hear the person with the microphone transmitter but also sounds closer to you. You could have the microphone transmitter picking up the television but still hear someone speaking to you nearby through the receiver microphone. You can easily adjust the sound from either microphone through volume controls on the receiver. The receiver's display shows the volume level, battery status and the selected channel.

- Built-in omni- and uni-directional microphones in both the transmitter and the receiver

- Sound is transmitted digitally
- Coded transmission for better confidentiality
- Choice of digital channels
- Supplied with a television connection lead
- Built-in rechargeable batteries with up to 20 hours' use from four hours' charge
- Operating range of up to 25m
- Size: 100 x 45 x 25mm

Options available

L320 Contego with neckloop

£435.00 exc VAT

£522.00 inc VAT

L321 Contego with headphones and earbuds

£435.00 exc VAT

£522.00 inc VAT

Accessories

L320A Optional neckloop

£46.51 exc VAT

£55.81 inc VAT

L334A Optional stethoset

£14.48 not VAT exempt

L321B Optional earbuds

£8.91 not VAT exempt

L320C Additional receiver

£312.35 exc VAT

£374.82 inc VAT



Watch this face

Lipreading is a valuable communication skill that can help you feel less isolated. *Watch this face* contains photographs of lipshapes that illustrate sounds, exercises to practice and common pitfalls so that you can avoid mistakes. An invaluable back-up to lipreading classes, the book can be used alone with a mirror, with a friend or with the TV. It will give you all the help, reassurance and confidence you need to start communicating better.

RN/3585/04/10

Watch this face

£8.99 not VAT chargeable

Introduction to the DigiSystem

This is our top-of-the-range system. The key difference between the DigiSystem and all of our other listeners is that you can pick up sound from several microphones. With a standard listener you just have the receiver microphone. With our other elite listeners you can pick up sound from the transmitter and receiver microphones, but with the DigiSystem you can

listen to more than one transmitter microphone. This means it is easier to listen to several people, making it particularly useful in meetings, classrooms, courtrooms, or wherever several people are talking.

The DigiSystem is extremely flexible because you can mix and match microphones. The sound quality is exceptional

because it can distinguish between voices and other noise, so it removes background sounds and enhances speech. A complete DigiSystem must consist of the DH-10 receiver and at least one microphone.

The DigiSystem is a high spec technical product so please contact our Customer Service



DigiSystem receiver

DH-10 Handset receiver

This discreet receiver comes with a neckloop for hearing aid wearers with a T setting and earbuds for people without hearing aids.

- Volume control and stereo sound
- You can set the volume separately for each ear
- Built-in rechargeable battery with 10 hours use from two hours charge
- Mains charger included
- Operating range of up to 30m
- Size: 74 x 39 x 13mm

L364 DH-10 Handset receiver

£695.00 exc VAT

£834.00 inc VAT



DM-20



DigiSystem master transmitter microphones

These microphones control the sound from other microphones and transmit to one or more receivers. Without one of these you cannot listen to more than one microphone. They can also be used as a single microphone if desired.

DM-20 Personal microphone

An advanced lecture microphone with a built-in receiver which easily attaches to the clothing of the person talking.

- Push to talk function prevents the receiver from hearing any other microphone
- Smaller external microphone provided in pack
- Mute button that can prevent the receiver from hearing this microphone

- Automatic gain control reduces the impact of loud sudden noises, as well as increasing the volume of quieter sounds
- Built-in rechargeable battery with seven hours' use from two hours' charge
- Size: 84 x 39 x 14mm

L381 DM-20 Personal microphone

£1,105.00 exc VAT
£1,326.00 inc VAT

DC-20 Conference and Bluetooth mobile microphone

This microphone has all the features of the DM-20 plus three advanced options:

- Tone control equaliser

- It can also be connected to mobile telephones via Bluetooth and to landline telephones via the Telelink accessory (purchased separately - see right)
- The microphone has four options: Omni-wide (for listening in all directions), Omni-narrow (for listening to the front), Zoom (for more targeted listening) and Music (enhances musical notes)
- Stand supplied
- Size: 115 x 62 x 14mm

L365 DC-20 Microphone

£2,705.00 exc VAT
£3,246.00 inc VAT

Suggested packs:

Meeting room pack

Ideal for presentations whilst hearing questions from the rest of the room. It can be expanded with extra DM-10s.

- 1xDH-10 Receiver
- 1xDM-20 Personal microphone
- 1xDM-10 Personal microphone

L407 Meeting room pack

£2,495.00 exc VAT **£2,994.00** inc VAT

Conference and mobile phone pack

Hear more clearly during meetings and presentations, Bluetooth microphone can connect to your Bluetooth mobile phone.

- 1xDH-10 Receiver
- 1xDC-20 Conference microphone
- 1xDM-05 Personal mini microphone

L409 Conference and mobile phone pack

£3,890.00 exc VAT **£4,668.00** inc VAT



DM-10



DigiSystem slave transmitter microphones

These slave microphones can be used with the L364 DH-10 receiver (see page 45) on their own or they can be used as additional microphones if you use one of the master microphones.

DM-10 Personal microphone

This has the same specification as the DM-20 microphone except it is not a master microphone.

- Size: 84 x 39 x 14mm

L380 DM-10 Personal microphone

£695.00 exc VAT
£834.00 inc VAT



DM-05

DC-10 Conference and Bluetooth microphone

This has the same specification as the DC-20 microphone except it is not a master microphone.

L372 DC-10 Conference and Bluetooth microphone

£2,005.00 exc VAT
£2,406.00 inc VAT

Accessories

Telelink

You can listen to your telephone hands-free with this sleek microphone transmitter. Can only be used with DC-20.

L367 Telelink

£371.00 exc VAT
£445.20 inc VAT

DM-05 Personal mini microphone

This is one of the world's smallest microphones for digital sound transmission. It has a similar specification to the DM-10 (left), although there is no mute button. It can be clipped onto clothes or worn around the neck.

- Size: 56 x 26 x 9mm

L379 DM-05 Personal mini microphone

£490.00 exc VAT
£588.00 inc VAT

Conference and mobile/landline phone pack

This pack will help in many common situations where you'll need hearing assistance.

- 1xDH-10 Receiver
- 1xDC-20 Personal microphone
- 1xDM-05 Personal mini microphone
- 1xTeleLink

L408 Conference and mobile/landline phone pack

£4,261.00 exc VAT **£5,113.20** inc VAT



Multi-purpose neckloop

This neckloop normally accompanies our Sonido (see page 36). You can also use it with many other products, including MP3 players, hi-fis and radios, by using the 3.5mm jack socket. The neckloop has a 0.9m lead and transmits sound to your hearing aid via the T setting.

L335A Adult neckloop

£12.50 exc VAT

£15.00 inc VAT



iPod not included

iLoop

This simple and modern hands-free neckloop has been designed for use with MP3 players, Walkmans and CD players. Simply place it around your neck and plug in. It is the perfect partner when you are on the move. The lightweight, battery-powered unit provides an on/off indicator light switch with volume control, offering

clearer sound quality for hearing aid wearers with a T setting.

- Cable length of 1m with 3.5mm jack plug
- Green LED when switched on

L309 iLoop

£32.50 exc VAT

£39.00 inc VAT

2 x AAA batteries included





Silent headphones for hearing aid wearers only

This model offers 'silent' inductive loop pick-up for use with a hearing aid only.

- Supplied with 1.15m lead and 3.5mm stereo jack plug and adjustable headband

M416 Silent headphones

£21.95 exc VAT

£26.34 inc VAT



Ear hooks for easier listening

Enjoy clearer sound from personal listeners, MP3, DVD or CD players with these lightweight inductive ear hooks. Simply place the hooks behind your hearing aid and they will transfer sound directly to the aid when on the T setting. You enjoy clear sound because they are close to your ear. Unlike with other earhooks you can also use these with your mobile phone. A single ear hook with an attached ear bud is perfect for single hearing aid wearers, who can now listen to music in both ears.

- Approx. 1m lead length
- Supplied with 2.5mm jack plug with 3.5mm adaptor
- Works with a variety of audio devices
- Powered by your audio device
- **An adaptor for your mobile may be required - please call Customer Service**

01733 361199 (telephone) or

01733 238020 (textphone)



Options available

M479 Mobile and audio single ear hook

£24.99 exc VAT

£29.99 inc VAT

M480 Mobile and audio dual ear hooks

£24.99 exc VAT

£29.99 inc VAT

M483 Mobile and audio ear hook and bud

£24.99 exc VAT

£29.99 inc VAT



Megaloop

The Megaloop is a powerful and flexible loop system. There are two microphone sockets. One is used when the microphone or scart connection is fixed to a TV or hi-fi speaker and the other can help with environmental sounds like the phone or doorbell or for conversation. We recommend that the cable is laid around the perimeter of the room at skirting board height.

- Supplied with 38m of loop cable
- Supplied with one microphone and scart lead
- Volume and tone control
- Headphone socket
- Size: 50 x 186 x 82mm

L295 Megaloop

£82.00 exc VAT

£98.40 inc VAT

Accessories

L295B Megaloop pad with 10m cable

£25.00 exc VAT

£30.00 inc VAT

L295C Additional microphone with 1.5m lead

£18.00 not VAT exempt

L295D Replacement loop cable

£18.00 not VAT exempt



Megaloop+

Developed from the Megaloop, this loop comes with two built-in scart sockets, microphone and line inputs. This means you can connect it to your video, satellite system and DVD in a number of possible ways without any extra junction boxes. You can even use it when the TV scart sockets are already in use. A scart fix feature is also available which solves the problem of hearing the 'wrong' sound by a simple press of a button. This system is easy to install and comes with all the leads you need.

- Supplied with 38m of loop wire and one microphone
- Volume and tone control
- Two microphone sockets
- Size: 70 x 181 x 83mm

L301 Megaloop+

£99.00 exc VAT

£118.80 inc VAT

Prompt delivery. Easy to install, clear instructions. Amplifies television sound well.

Mr and Mrs Ward, Colchester



EchoSL

This product is the only domestic loop system with a remote control. The remote controls the volume, tone, mute function and power.

If you use your hearing aid's T setting when using your telephone, then this product has an extra benefit. You can simply mute the EchoSL with the remote rather than getting up to switch off the loop manually (as you would have to with other models) in order to avoid interference.

The EchoSL comes supplied with a stand so that it can be stood upright or flat, to suit the space available around your TV. The remote won't work with other loops, but the EchoSL will work with your existing loop wire. This means you can simply swap this amplifier with your current one without removing and refitting the wire.

- Supplied with 38m of cable
- Includes a remote control
- Scart and phono leads supplied
- Size: 45 x 190 x 81mm

L406 EchoSL

£125.00 exc VAT

£150.00 inc VAT

2 x AA batteries included



LA215 Loop amplifier

This loop system has a modern, flexible design. You can stand the unit upright or lay it flat. It is ideal for listening to the television, radio or hi-fi and it includes a microphone for sound pick-up and an audio connection lead kit for direct connection to the TV.

- Supplied with 37m of loop cable
- Adjustable volume and tone controls
- Scart and phono-phono leads supplied
- A switch to select which device you are listening to, such as the stereo or TV
- Size: 40 x 180 x 150mm

L337 LA215 Loop amplifier

£81.50 exc VAT

£97.80 inc VAT

Accessories

X672 Replacement tie-clip microphone with 4m lead

£11.88 not VAT exempt

X675 Optional loop pad

£19.50 exc VAT

£23.40 inc VAT

X702 3m Extension cable for microphone

£5.40 not VAT exempt





LoopHear

This is an easy-to-use rechargeable, portable loop system that enables people with hearing aids with a T setting to hear conversations more clearly. Because of its metre range, LoopHear is ideal for one-to-one applications such as at service counters, railway stations, banks and in small meetings. The lightweight unit incorporates an induction loop and a built-in microphone and also has a useful base plate which can be fixed to a surface for a more permanent set-up.

- Socket to add an extra microphone
- Volume control
- NiMH battery is rechargeable and replaceable
- Removable base plate for a semi-fixed position
- Outer carton of packaging can be used as a handy carrying case
- Power on LED
- Low battery indication (flashes red)

- Battery time of eight hours operation time or two days standby time
- Battery power lasts for four months when unit is switched off
- Operating range of up to 1m
- Supplied with stethoset and lapel microphone
- Headphone socket
- Size: 210 x 240 x 70mm

L371 LoopHear

£112.00 exc VAT

£134.40 inc VAT



Ezee Loop

One of our most popular products, the Ezee Loop is an extremely flexible, portable system. It is perfect for one-to-one discussions, counters or reception desks. The Ezee Loop's built-in microphone picks up the speaker's voice and relays it to people's hearing aids. It has a range of about 1m. Using the optional boundary microphone it can pick up sound at a greater distance, such as in a small meeting room. It can be connected to a television or video conferencing equipment with a scart lead.

- Powered by a rechargeable internal battery or mains (battery included)
- The battery gives a minimum of 24 hours' use
- Operating range of up to 1m
- Size: 250 x 280 x 110mm

PL13 Ezee Loop

£115.00 exc VAT

£138.00 inc VAT

Accessories

PL13A External boundary microphone lead 2m

£35.81 not VAT exempt



Portable loop system

These portable loops can provide hearing support for whole rooms – up to 160m². They reduce background noise, making it easier for people who are deaf or hard of hearing to hear through their hearing aids. The system comprises at least one microphone, an amplifier and a length of cable that is looped around the room. Although we recommend installed loops in many cases, portable systems can be useful as the lower cost and flexibility make them ideal for smaller organisations or people on the move.

- Size: 130 x 460 x 330mm

Prices range from around **£1,000** (exc VAT) to over **£1,800**

Contact us for further details.



Portable infrared system

Portable infrared systems are an increasingly popular alternative to loop systems. Unlike a loop, there is no need to run a cable around the room, making it quicker to set up. An invisible infrared signal is sent to a receiver worn by the user. This means infrared is the best option when confidentiality is an issue and where there is more than one system being used nearby, as no-one outside the room can listen to the conversation. You do not have to be a hearing aid wearer to benefit from an infrared system.

- Size: 270 x 500 x 540mm

Our portable infrared systems start around **£2,050** (exc VAT) and can cover an area up to 220m² (for the largest system).

Contact us for further details.



Clipboard

This product enables organisations to provide hearing support for hearing aid wearers. The clipboard is rechargeable, portable and lightweight, making it ideal for home visits, sales assistants in shops, or when taking orders in a restaurant. The clipboard contains an internal induction loop and an integral microphone, which picks up the voice of the user and relays it directly to a hearing aid wearer using the T setting.

- No installation required
- On/off switch and visual indicator showing power, signal strength and battery life
- Automatic standby
- Built-in rechargeable battery and charger supplied
- Battery life of up to seven hours
- Operating range of up to 1.2m
- Size: 320 x 230 x 30mm

L312 Clipboard

£144.00 exc VAT

£172.80 inc VAT