



## User Manual, Bed-shaker (GB)

---

*Bellman Bed-shaker can be used with all the Bellman Visit receivers including the Pager when mounted in the charger, as a means of waking up with help of vibrations.*

### This is what to do!

1. Connect the Bed-shaker (1) by inserting the connector (2) in the indicated output (3) on the receiver (see above text).
2. Press the test button on the transmitter you have for the Visit System and check that the bed vibrator gives a signal.

Note! When you connect the bed-shaker to the Bellman Visit Pager must be inserted in the charger for the bed-shaker to function.

### Production description

- 1) Slip protection
- 2) Connector

